

Anti-Human CD8-Phycoerythrin

Product Information

Part Number: 25168-00

Size: 100 tests

Volume per Test: 10µl

Clone: RPA-T8

Species reactivity: Human

Isotype: Mouse IgG1

Storage Buffer: Aqueous buffered solution containing 4.6% BSA, 5% Trehalose and 0.01% Sodium Azide

Absorption Maximum: 565nm

Emission Maximum: 575nm

Preparation and Storage

Antibody was purified by Protein-Afinity chromatography

Storage temperature: 4°C and protected from prolonged exposure to light

Shelf Life: **One year from date of receipt**

Application Notes

This product has been reported to work in the following applications. This information is derived from testing within our laboratories, and peer-reviewed publications.

- Flow Cytometry
- Immunohistochemistry (frozen)

ReaMetrix India Pvt. Ltd.

www.reametrix.com

For further information and pricing, please call us at:

India

+91 80 2837 8693/5

United States

+1 650 620 9253